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## Experiment 1 Extraction, Isolation and Detection of Anthraquinones from *Rheum palmatum* L.

### 1. Objective

1. To master extraction, isolation principles and methods of anthraquinones from *Rheum palmatum* L.
2. To master the separation principle and method of silica gel column chromatography, and to learn dry column-packing method.

3. To learn the method of preparation of silica gel thin layer plate.
4. To learn the chemical methods of detecting anthraquinones.

## 2. Experimental principles

*Rhubarb* was recorded in *Shen Nong's Herbal Classic*. It has the effects of reducing internal heat, defaecation, cooling blood, detoxification, eliminating stasis, and dredging collateral channels. It is traditionally used for defaecation, promoting digestion, clearing heat, and detoxication as a folk medicine. *Rhubarb* is a famous plant purgating drug, and has been recorded early in pharmacopoeia of many countrys. All the roots and rhizomes of *Rheum palmatum* L., *Rheum officinale* Baill, and *Rheum tanguticum* Maxim. et Regll could be used as *Rhubarb*. Its chemical constituents contain various hydroxy-anthraquinones and their glycosides. The main anthraquinone aglycons are listed in the following table:

Table 10-1 Characteristics of main anthraquinone aglycons

R <sub>1</sub>	R <sub>2</sub>	Compounds	Crystal form	Melt point(°C)
-H	-COOH	Rhein	Yellow needle crystal	318~320
-CH <sub>3</sub>	-OH	Emodin	Orange needle crystal	256~257
-H	-CH <sub>2</sub> OH	Aloe-emodin	Orange needle crystal	206~208
-CH <sub>3</sub>	-OCH <sub>3</sub>	Physcion	Brick-red needle crystal	207
-H	-CH <sub>3</sub>	Chrysophanol	Golden flaky crystal	196

Polarity of anthraquinone aglycons in *Rhubarb* is different because of their different structures, especially substitution groups. Their polarity, from strong to weak, is in the sequence of Rhein, Emodin, Aloe-emodin, Physcion and Chrysophanol.

The method of acid hydrolysis of *Rhubarb* Powder is used for higher aglycons' yield. The hydrolysis product, aglycon, is extracted with ethyl ether, and the ethyl ether extraction is separated by silica gel column chromatography according to their polarity.

## 3. Instruments and reagents

### 1. Experimental instruments

Glass chromatographic column (21 mm×200 mm, 24#), globe-shaped funnel (250 ml, 24#), flask (100 ml), evaporating dish, test tube, glass thin layer plate(5×20 cm), filter paper, glassfunnel, iron support stand, absorbent cotton, round bottom flask, economical allihn condenser, thermostat water bath, Büchner flask, Buser funnel, ultraviolet lamp.

### 2. Reagents and drugs

*Rhubarb* Powder, silica gel (200-300 mesh), silica gel G, ethyl ether (AR), petrolum ether (60-90°C, AR), ethyl acetate (AR), 20% H<sub>2</sub>SO<sub>4</sub>, 0.5% CMC, 1% alizarin

alcohol solution, 0.5% magnesium acetate methanol solution.

#### 4. Experimental procedures

##### 1. Extraction and Isolation of Anthraquinone from *Rheum palmatum* L.

The dried and finely ground plant material (5 g) of *Rheum palmatum* L. is hydrolyzed with 20% H<sub>2</sub>SO<sub>4</sub> (60 ml) at water bath for 2 hours (operating in 250 ml round bottom flask), hydrolyzed material is obtained after vacuum filtration (washed to neutral with water) and air drying.

Dried hydrolyzed material is extracted with ethyl ether (25 ml) at room temperature about 2-3 days. The ethyl ether extraction is subjected to separation over a silica gel column and eluted with a petroleum ether-EtOAc gradient solvent system (7:3, ~200 ml). Eluents are collected with test tubes, 5~7 ml per fraction. The fractions are defined by the results of TLC (silica gel, developing solvent: petroleum ether-EtOAc 7:3). Each fraction then is evaporated and obtains pure compounds. Column chromatography operation is as follows:

The ethyl ether extraction is adsorbed onto 0.5~1.0 g silica gel (200~300 mesh) in evaporating dish and the solvent is evaporated at water bath, thus obtained sample silica gel. Dry column-packing method: silica gel (30 g, 200~300 mesh) is added to the chromatography column, then sample silica gel is packed on the top, some silica gel and absorbent cotton are placed on the top of the sample to prevent washing the sample in a mess when adding solvent on the top.

##### 2. Chemistry detection reactions

###### (1) Bornträger's reaction

Sample eluent and 1% alizarin alcohol solution each 2 ml are added 2 ml 1% NaOH at test tube, respectively, shake and then standing, observe color changes.

###### (2) Magnesium acetate test

Taking a drop of sample eluent and 1% alizarin alcohol solution are added on the same filter paper strip (3×10 cm), and then 0.5% magnesium acetate methanol solution is sprayed on this strip, observe color changes after heated at 80°C about 5 min.

#### 5. Questions

1. What are the principles of extraction and isolation of anthraquinones from *Rheum palmatum* L.?
2. Try to write extraction and isolation flow chart of anthraquinone aglycons.
3. What about the R<sub>f</sub> values of five anthraquinone aglycons when they are separated by silica gel thin-layer chromatography and developed with petroleum ether: ethyl acetate (7:3), and why?

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## Experiment 3 Determination of Rutin and Quercetin by UV Spectroscopic Method

### 1. Objective

- (1) To master the relationship between UV spectra and the structures of flavonoids.
- (2) To determine the structures of flavonoids by the UV spectra.
- (3) To master the operation method of T6-ultraviolet spectrometer.

## 2. Experiment principles

There are two main absorption peaks, called bands I and II, respectively, for the majority of flavonoids in methanol scanning from 200 to 400 nm in ultraviolet band, caused by two cross-shaped conjugation systems i. e. cinnamyl and benzoyl groups (Figure 10-1). The positions and shapes of the bands I and II are influenced by substituents on rings A and B. So according to the position, shape or intensity of bands I and II as well as the transformation of positions of the band with some diagnostic reagents, we can investigate the structures of flavonoids.

Common diagnostic reagents include sodium methoxide, sodium acetate, sodium acetate-boric acid, aluminium trichloride and aluminium trichloride-hydrochloric acid.


Figure 10-1 The structures of benzoyl and cinnamyl groups in flavonoids

Procedures of UV spectral determination in this experiment:

- a. Determination of the sample's UV spectrum in methanol.
- b. Determination of the sample's UV spectra after adding every diagnostic reagent in methanol solution, respectively.
- c. Deduce the structure of the sample according to the above spectra.

## 3. Instruments and reagents

- (1) Experiment instruments: T6-Ultraviolet spectrometer.
- (2) Reagents and drugs

Rutin sample, weigh precisely 10 mg of rutin into a measuring flask and dissolve with methanol to 100 ml. Pipet above solution 5 ml into a 50 ml measuring flask and dilute with methanol to 100 ml exactly, e.i. 10  $\mu\text{g/ml}$ .

Quercetin sample, weigh precisely 10 mg quercetin to a measuring flask and dissolve with methanol to 100 ml. Pipet above solution 5 ml into a 50 ml measuring flask and dilute with methanol to the volume accurately, 10  $\mu\text{g/ml}$ .

Sodium methoxide, weigh sodium metal 2.5 g and cut it into small pieces, then

gradually add the pieces into 100 ml methanol. The prepared solution must be stored in the airtight glass container and standby application.

Anhydrous sodium acetate, add sodium acetate in evaporating dish and dry at 120 degrees for 2 hours, then grind and standby application.

Anhydrous boric acid, grind and standby application.

Anhydrous aluminium trichloride, weigh about 1 g anhydrous aluminium trichloride and dissolve with 20 ml chromatogram class methanol, stored for 24 hours.

Hydrochloric acid, take 50 ml concentrated hydrochloric acid and diluted with 100 ml distilled water, blending and standby application.

#### 4. Determination procedures

(1) Determination of rutin UV spectra (scanning from 220 to 550 nm)

1) Determination of rutin UV spectra in MeOH.

2) Determination of rutin UV spectra in MeOH + NaOCH<sub>3</sub>: Load test sample to silici-cuvette, add 3 drops sodium methanol solution to the silici-cuvette and then scan immediately.

3) Determination of rutin UV spectra in MeOH + AlCl<sub>3</sub>: Load test sample to silici-cuvette and add 6 drops AlCl<sub>3</sub> solution to the silici-cuvette and then scan.

4) Determination of rutin UV spectra in MeOH + AlCl<sub>3</sub> + HCl: Adding 3 drops of hydrochloric acid solution to the above sample solution, and then scan immediately.

5) Determination of rutin UV spectra in MeOH + NaOAc: Load test sample to silici-cuvette and add anhydrous sodium acetate, shake till the sediment on the bottom reached ca. 2 mm, and then scan.

6) Determination of rutin UV spectra in MeOH + NaOAc + H<sub>3</sub>BO<sub>3</sub>: Adding enough hydrochloric acid solution to the above sample solution, and then scan immediately.

(2) Determination of quercetin UV spectra (scanning from 220 to 550 nm):

All the determination methods refer to rutin as described above.

#### 5. Data analysis

(1) To explain the relationship between UV spectral data and the structure of rutin.

(2) To explain the relationship between UV spectral data and the structure of quercetin.

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